



KACEY DIAGNOSTICS's

ORGANO-PHOSPHATE Diagnostic Poison Test Strips

“The only ONE step, SIX minute test for Canine and Feline use!”


KCOP diagnostic strips detects the presence of Organophosphates commonly found in over 175 insecticides, pesticides, cholinesterase inhibitors and OP Carbamate products. The only veterinary test strip to provide STAT, in house, quantitative measurement of OP in plasma for both Canine and Feline in six minutes. When time counts, KCOP strips are ready to make the difference.

In House, KACEY “OP” TEST STRIPS	OFF site method
1 drop of plasma (20ul) pulled in clinic	Blood drawn and sent out
1 step process with results in 6 minutes	Earliest results in 24 hours
Quick color chart comparison	New tests must repeat process

Specifications:

- Feline sensitivity level 75 U/dl @ NLT
- Canine sensitivity level 125 U/dl @ NLT
- Only 20ul of plasma needed
- Results in 6 minutes
- Easy color comparison chart
- 14 months expiration
- Low unit cost per strip

KCEGT: Simplicity and Speed!



KACEY “OP” strips are packaged in either 5 strip or 10 strip bottles for convenience. Easy to use 20 ul pipettor and pipettes are also available. Contact us for immediate service or the name of your nearest distributor.

Contact us or your local Kacey Diagnostics dealer for more information about the KCOP test strips.

Authorized Dealer:



KACEY™
DIAGNOSTICS

1854 A Hendersonville Rd
Asheville, NC 28803

Ph: 1(828)685-3569

Fax: 1(828)685-3571

www.KaceyDiagnostics.com

“Kacey Products Providing Better Animal Care”